AUDIT® MicroDROPTM Urine Dipstick Control

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ENGLISH

INTENDED USE

Audit® MicroDROPTM Urine Dipstick Control is a reference control consisting of human source materials fortified to target levels with compounds that produce the desired reaction when tested by the methods indicated. It is intended to simulate human patient urine samples for the purpose of monitoring the precision of the Multistix®, Uristix®, Chemstrip®, and vChem™ urinalysis reagent strips, the Clintek Status® urine chemistry analyzers, as well as for confirmatory tests such as Acetest® and Clinitest®.

SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Audit® MicroDROPTM Urine Dipstick Control may be used as one would use human urine to obtain the stated values.

REAGENTS

Audit® MicroDROPTM Urine Dipstick Control is an in vitro diagnostic control composed of a buffered bacteriostatic and fungistatic human serum matrix. Audit® MicroDROPTM Urine Dipstick Control is a liquid product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBSAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

Audit® MicroDROPTM Urine Dipstick Control is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

STORAGE AND STABILITY

Audit® MicroDROPTM Urine Dipstick Control is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for urinalysis reagent strips, Audit® MicroDROPTM Urine Dipstick Control has an open vial stability of up to 30 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

MATERIALS PROVIDED

Audit® MicroDROPTM Urine Dipstick Control, 10 x 5 mL

PROCEDURE

Follow the manufacturer's instructions provided for urinalysis reagent strips. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

INSTRUCTIONS FOR USE

1. The Audit® MicroDROPTM Urine Dipstick Control should be stored at 2-8°C. Do not shake. Do not mix mechanically.
2. Remove the controls from the refrigerator and allow to come to room temperature (20-25°C) at least 15 minutes, depending on the remaining volume. Mix gently by inversion to assure complete mixing of the contents.
3. Remove cap and invert bottle. While holding the dipstick, gently squeeze the dropper bottle, and touch the tip of the bottle to the dipstick, thoroughly saturating each pad. Turn dipstick on its side to drain excess control material onto an absorbent material.
4. Read the Urine Dipstick in accordance with the manufacturer's instruction manual.
5. After sampling, wipe off dropper tip, replace cap and return to 2-8°C when not in use.

EXPECTED VALUES

The expected ranges for visual readings for each level, based on interlaboratory data established by comparing the dipstick reaction that occurs with the controls to the color comparison chart with multiple lots of each manufacturers' dipsticks or reagent tablets, is

provided below.

LIMITATIONS OF THE PROCEDURE

Some analytes may take longer to develop than the indicated times. Bilirubin and urobilinogen may not show characteristic colors on some strips. Make sure that the vial is brought to room temperature before testing. If the liquid in the vial becomes frozen, discard and use another vial, as the results will not be valid.

Dispose any discarded materials in accordance with the requirements of your local waste management authorities.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K064M-10	Urine Dipstick Control	10 x 5 mL

Distributed by AUDIT MicroControls, Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

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