AUDIT[®] MicroCV[™] Procalcitonin Method Validation Set

REF K842M-25

25 x 1 mL

LOT 06289

Aalto Scientific Ltd
Carlsbad, CA 92008
USA

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ENGLISH

INTENDED USE

The Audit[®] MicroCV[™] Procalcitonin Method Validation Set is assayed quality control material consisting of twenty-five levels of bovine based serum albumin. Levels 1 - 22 contain the following analyte: Procalcitonin. Three vials of negative samples are provided. It is intended to simulate human patient serum samples for purpose of determining calibration verification and verification of reportable range for Procalcitonin.

This product is intended for use with quantitative assays on the indicated analyzer provided in the labeling. The Audit[®] MicroCV[™] Procalcitonin Method Validation Set should not be used for calibration or standardization of the Procalcitonin assays. The Audit[®] MicroCV[™] Procalcitonin Method Validation Set is for In Vitro Diagnostic use only.

SUMMARY AND PRINCIPLE

As defined in the Clinical Laboratory Improvement Amendments of 1988 (CLIA) by the Centers for Medicare and Medicaid Services (CMS) and the Centers for Disease Control (CDC), each laboratory must demonstrate that it can obtain performance specifications for an FDA-cleared or approved test system comparable to those established by the manufacturer for accuracy, precision and reportable range of test results for the test system¹. Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Audit[®] MicroCV[™] Procalcitonin Method Validation Set may be used as one would use human serum to verify these performance specifications.

WARNINGS AND PRECAUTIONS

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Audit[®] MicroCV[™] Procalcitonin Method Validation Set is intended solely for the purpose of in vitro diagnostic use as described on the label. Audit[®] MicroControls[™], Inc. will not be liable for any unclaimed damages arising from any other usage.

MATERIALS PROVIDED

Audit[®] MicroCV[™] Procalcitonin Method Validation Set, 25 x 1 mL

STORAGE AND STABILITY

Audit[®] MicroCV[™] Procalcitonin Method Validation Set is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for its assays, Audit[®] MicroCV[™] Procalcitonin Method Validation Set has a reconstituted stability of up to 5 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

PROCEDURE

Follow the manufacturer's instructions provided for establishing and verifying performance specifications for accuracy, precision and reportable range of test results for the test system. Verify that the lot number on each vial matches the package insert. To avoid evaporation, do not leave the vial uncapped. Q.C. requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements. Method Validation material should be run:

1. before the introduction of a new test system into routine use
2. whenever conditions change for which a method has been validated, e.g., an instrument change with different characteristics
3. whenever the method is changed and the change is outside the original scope of the method.

INSTRUCTIONS FOR USE

1. Remove a vial from the package, reconstitute with 1 mL of deionized water using a volumetric pipet, and gently swirl occasionally for 10 minutes. Do not shake. Do not mix mechanically. The Negative vials are liquid and do not need to be reconstituted.
2. Refer to instrument or assay instruction manual for establishing and verifying performance specifications for accuracy, precision and reportable range of test results for the test system.
3. After sampling, replace stopper and return to original package at 2-8°C to obtain the maximum open vial stability.

CALCULATIONS OF RESULTS

Refer to the manufacturer's instructions and guidelines for reporting data for the analysis of the test system's method validation studies.

LIMITATIONS OF THE PROCEDURE

Make sure that each vial is brought to room temperature before testing. If the contents of any of the vials become frozen, discard all vials and request a replacement set, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

EXPECTED VALUES

The analyte concentrations in this insert were derived from multiple replicate analyses. Actual results obtained may vary depending on instrumentation, methodology and assay temperature. Results may also be dependent on the accuracy of the instrument/reagent system calibration. The degree of acceptable non-linearity is an individual judgment based on methodology, clinical significance and medical decision levels of the test analyte. The material and information presented here in no manner constitutes an overruling of any federal, state or other regulatory body's regulations and/or guidelines.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K842M-25	Procalcitonin Method Validation Set	25 x 1 mL

Distributed by AUDIT MicroControls, Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

¹Federal Register 42 CFR Part 493, Department of Health and Human Services, January 24, 2003; §493.1253, (b)(1)(i)(A), (b)(1)(i)(B) and (b)(1)(i)(C).

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