



Audit® MicroCV™ Tumor Markers Linearity Set For Use With Siemens Centaur® Immunoassay Systems

Cat. No. **K805M-5**
Contents 10 x 1mL

Lot No. 06195A, 06195B, 06195C,
06195D, 06195E
06195AA, 06195BB, 06195CC,
06195DD, 06195EE

For In Vitro Diagnostic Use Only.

©Audit® MicroControls™, Inc., Las Vegas, NV 89102, (866) 252-8348

Expires 05/03/2012

INTENDED USE

Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems consists of five levels of human based serum. Set #1 contains the following analytes: Alpha fetoprotein (AFP), Carcinoembryonic antigen (CEA), Prostate-specific antigen (PSA), and CA 125. Set #2 contains CA15-3. These five levels of each set (#1 and #2) demonstrate a linear relationship to each other for their respective analytes, reagents and instruments¹.

This product may also be used as unassayed quality control material for these analytes. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

SUMMARY AND PRINCIPLE

As defined in the Clinical Laboratory Improvement Amendments of 1988 (CLIA) by the Centers for Medicare and Medicaid Services (CMS) and the Centers for Disease Control (CDC), each laboratory must revalidate each test method's analytical measurement range (AMR) at least every six months as well as following changes in lots of analytically critical reagents or major system components². Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems may be used as one would use human serum to verify and validate the AMR.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBSAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit® MicroControls™, Inc. shall not be liable for any unclaimed damages arising from any other usage.

STORAGE AND STABILITY

Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be reconstituted immediately and used according to the instrument manufacturer's instructions.

When used for the purpose of verifying and validating the AMR for non-waived tumor marker testing methods, it is recommended that Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems be used within seven days after reconstitution (four days for CA 15-3) when stored tightly capped at 2-8°C.

Leaving the vial uncapped, or prolonging its time at room temperature, will void this reconstituted stability claim. Make sure the contents of the vial are well mixed before use.

PROCEDURE

Follow the manufacturer's instructions provided for quality control and for verifying and validating the AMR. Verify that the lot number on each vial matches the package insert. To avoid evaporation, do not leave the vial uncapped. Q.C. requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements. CLIA directed calibration verification linearity material should be run³:

1. every six (6) months.
2. when a complete change of reagents for a procedure is introduced.
3. there is major preventive maintenance or replacement of critical parts that may influence test performance.
4. when control materials reflect an unusual trend or shift, or are outside of the laboratory's acceptable limits.
5. when the laboratory's established schedule for verifying the reportable range for patient test results requires more frequent calibration verification.

Materials provided

- Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems, 10 x 1mL

INSTRUCTIONS FOR USE

1. Remove a vial from the package, reconstitute with 1ml of deionized water, and gently swirl occasionally for 10 minutes. Do not shake. Do not mix mechanically.
2. Refer to instrument or assay instruction manual for verifying and validating the AMR.
3. After sampling, replace stopper and return to original package at 2-8°C to obtain the maximum seven-day reconstituted stability.

CALCULATIONS OF RESULTS

Each set of Audit® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems is prepared in a manner such that an equal distance exists between each consecutive level. This dilution scheme is consistent with the NCCLS recommendation¹ for preparing linearity sets.

Once each vial of the total set is tested, raw data may be entered via the AUDITOR™ QC Program at www.auditmicro.com. An on-line graph showing actual values versus predicted values for each analyte is then available to print, along with slope and intercept data. Call (866) 25-AUDIT for more information.

LIMITATIONS OF THE PROCEDURE

Make sure that each vial is brought to room temperature before testing. If the contents of any of the vials become frozen, discard all vials and request a replacement set, as the results will not be valid.

¹Dilution schemes are based on guidelines provided by The National Committee for Clinical Laboratory Standards (NCCLS) in approved guideline EP6-A, "Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach, Approved Guideline", April 2003.

²Federal Register 42 CFR Part 493, Department of Health and Human Services, January 24, 2003; p.3691.

³Federal Register 42 CFR Part 493, Department of Health and Human Services, January 24, 2003; §493.1255, (b) (1) (ii).

LINEARITY FACTOR

Each lot of product is manufactured such that a linear relationship exists among levels. The analyte concentrations in this insert were derived from multiple replicate analyses. Actual results obtained may vary depending on instrumentation, methodology and assay temperature. Results may also be dependent on the accuracy of the instrument/reagent system calibration. The degree of acceptable non-linearity is an individual judgment based on methodology, clinical significance and medical decision levels of the test analyte. The material and information presented here in no manner constitutes an overruling of any federal, state or other regulatory body's regulations and/or guidelines.

Analytes include:

ORDERING INFORMATION

Product Number	Product Description	Product Packaging
K805M-5	AUDIT® MicroCV™ Tumor Markers Linearity Set for use with Siemens Centaur® Immunoassay Systems	10 x 1mL

Set 1 - Centaur®

	UNITS	A	B	C	D	E
AFP	ng/mL	5.80	229	458	658	910
CEA	ng/mL	1.49	24.1	52.0	73.1	91.4
PSA	ng/mL	0.11	25.6	52.8	89.1	119
CA125	U/mL	21.1	163	303	457	597

Set 1 – Centaur CP®

	UNITS	A	B	C	D	E
AFP	ng/mL	4.88	213	430	646	861
CEA	ng/mL	1.57	29.2	53.4	85.7	108
PSA	ng/mL	0.09	26.7	56.4	87.5	130
CA125	U/mL	22.2	172	312	458	603

Set 2 - Centaur®

	UNITS	A	B	C	D	E
CA15-3	U/mL	16.1	74.0	119	160	190

Set 2 – Centaur CP®

	UNITS	A	B	C	D	E
CA15-3	U/mL	13.7	69.9	117	160	194