



AUDIT® MicroControls™ recommends the following checklist if your initial calibration verification fails for any specific analyte or analytes:

1. Quality Control Material

- Are there patterns among controls (e.g. all are below the mean, all are above the mean)?
- Are there any noticeable trends or shift over time?
- How is the accuracy and precision?

2. Acceptable Range

- Re-examine your laboratory's determination of the acceptable range for calibration verification material.
- What is your laboratory's current range around the expected target value for that specific analyte in question?

3. Reagent Changes

- Have there been any changes to your reagent?
- New lot of reagent?
- Different manufacturer?
- New formulation of current reagent (check package insert)?

4. Instrument Maintenance Logs

- Review daily, weekly, monthly, quarterly, semi-annual and annual logs for any deviations or changes.

5. Environmental

- Has the instrument been moved recently?
- Any changes to the environment of the instrument and its surroundings?

6. Servicing

- Has the instrument been serviced recently?
- Any software or hardware upgrades or changes?

7. Operation

- Are there new instrument operators?
- Any recent modification to the technique in how the assay is run?

8. Comparative Method

- Is there another nearby laboratory that can also run the calibration verification material to compare results?
- If all the above has been performed, and there are still problems, re-calibrate the instrument.
- If the instrument still does not perform within laboratory control limits, call the instrument manufacturer for further troubleshooting.