

AUDIT MicroControls, Inc. Announces Audit® MicroCV™ Blood Gas Linearity Set

August 15, 2011 (CARLSBAD, CA), AUDIT MicroControls, Inc. is pleased to announce a new addition to our existing line of calibration verification/linearity products. The Audit® MicroCV™ Blood Gas Linearity Set has five unique levels of linearity material which allows for complete monitoring of that specific test system's reportable range, as defined by the manufacturer.

Packaged as 15 x 2 mL, this liquid has a shelf life of 1 year and an open vial stability of 60 seconds.

Recently adopted CLIA-88 regulations state that each laboratory revalidate each of its test method's analytical measurement range (AMR) at least every six months as well as following changes in lots of analytically critical reagents or major system components.

AUDIT MicroControls provides a free data reduction program to their customers, called AUDITOR™ QC. This service is a group of Internet based tools that perform statistical analysis calculations on the linearity data supplied by the lab and provides the end user peer group comparison graphs that may be printed or stored for future use.

Further information can be found at www.auditmicro.com.

About AUDIT MicroControls

Founded in 2002, AUDIT MicroControls is a wholly owned subsidiary of Aalto Scientific, Ltd. of Carlsbad, CA. Specializing in calibration verification/linearity products, and routine and esoteric chemistry quality controls, AUDIT products are available through a number of regional and national distributors. For more information about AUDIT, please call (866) 25-AUDIT or visit www.auditmicro.com.